

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT	MX 447 L	CONTROL L		(Exp.)		2024-08-05 (YYYY-MM-DD)
Rev 2						
LyseVet						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS Scil Vet ABC Plus Scil Vet ABC ESV 60				TOLERANCES TOLERANCE
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8				± 0.4
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.27				± 0.15
HB HGB	g/dl	5.9				± 0.4
	g/l	59				± 4
	mmol/l	3.66				± 0.25
HT HCT	%	16.3				± 2.0
	l/l	0.163				± 0.020
VGM MCV	µm ³ ; fl	72				± 4
TGMH MCH	pg	26.0				± 2.0
	fmol	1.61				± 0.12
CCMH MCHC	g/dl	36.1				± 3.0
	g/l	361				± 30
	mmol/l	22.42				± 1.86
IDR RDW *	%	15.0				± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	74				± 20
VPM MPV *	µm ³ ; fl	8.0				± 2.0
LYMPHO *	#	0.96				± 0.40
	%	53.5				± 8.0
MONO *	#	0.23				± 0.21
	%	13.0				± 11.7
GRANULO *	#	0.60				± 0.40
	%	33.5				± 10.0

* N/A for Micros Scil Vet ABC

LOT	MX 447 N	CONTROL N		(Exp.)		2024-08-05 (YYYY-MM-DD)
Rev 2						
LyseVet						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS Scil Vet ABC Plus Scil Vet ABC ESV 60				TOLERANCES TOLERANCE
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1				± 0.8
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.54				± 0.18
HB HGB	g/dl	13.2				± 0.5
	g/l	132				± 5
	mmol/l	8.20				± 0.31
HT HCT	%	36.8				± 2.5
	l/l	0.368				± 0.025
VGM MCV	µm ³ ; fl	81				± 4
TGMH MCH	pg	29.1				± 2.0
	fmol	1.81				± 0.12
CCMH MCHC	g/dl	35.9				± 3.0
	g/l	359				± 30
	mmol/l	22.29				± 1.86
IDR RDW *	%	14.5				± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	262				± 40
VPM MPV *	µm ³ ; fl	7.5				± 2.0
LYMPHO *	#	2.02				± 0.40
	%	28.5				± 5.0
MONO *	#	0.71				± 0.64
	%	10.0				± 9.0
GRANULO *	#	4.37				± 0.71
	%	61.5				± 10.0

* N/A for Micros Scil Vet ABC

LOT	MX 447 H	CONTROL H		(Exp.)		2024-08-05 (YYYY-MM-DD)
Rev 2						
LyseVet						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS Scil Vet ABC Plus Scil Vet ABC ESV 60				TOLERANCES TOLERANCE
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.7				± 1.6
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.73				± 0.20
HB HGB	g/dl	18.1				± 0.6
	g/l	181				± 6
	mmol/l	11.24				± 0.37
HT HCT	%	50.4				± 3.0
	l/l	0.504				± 0.030
VGM MCV	µm ³ ; fl	88				± 4
TGMH MCH	pg	31.6				± 2.0
	fmol	1.96				± 0.12
CCMH MCHC	g/dl	35.9				± 3.0
	g/l	359				± 30
	mmol/l	22.29				± 1.86
IDR RDW *	%	13.5				± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	498				± 45
VPM MPV *	µm ³ ; fl	7.4				± 2.0
LYMPHO *	#	2.17				± 1.00
	%	11.0				± 5.0
MONO *	#	1.77				± 1.59
	%	9.0				± 8.1
GRANULO *	#	15.76				± 1.97
	%	80.0				± 10.0

* N/A for Micros Scil Vet ABC