



Steroid Panel

LC-MS.

NEW



CE-IVDD solution for steroid determination in human serum and EDTA plasma by LC-MS.

THE NOVELTY IN TECAN'S REAGENTS

One single kit to quickly and reliably determine 18 different and clinically relevant steroids with only one extraction process.

BENEFIT WITH OUR STEROID PANEL LC-MS KIT*

LC-MS/MS is considered the gold standard in the determination of steroid hormones due to its high sensitivity!

- Steroid kit, CE-IVDD
- 18 steroids
- 2 chromatographic runs (guarantees the best possible performance of all analytes)
- Only one sample preparation with 96-well SPE-plate
- Manual or automated sample preparation possible (e.g. Tecan Resolvex® A200)
- All standards in one mix
- Individual isotopic standard for each analyte
- Total runtime (LC-MS): 10 min for 15 steroids, second injection for another 3 analytes in 6 min runtime

PRODUCTS

30191875

[Steroid Panel LC-MS >](#)

THE ANALYTES

11-deoxycorticosterone, 11-deoxycortisol, 17-hydroxypregnenolone, 17-hydroxyprogesterone, 21-deoxycortisol, aldosterone, androstenedione, corticosterone, cortisol, cortisone, dehydroepiandrosterone, dehydroepiandrosterone sulfate, dihydrotestosterone, estradiol, estrone, progesterone, testosterone, dexamethasone.



*Product and solution may not be available in all territories. Please consult your local TECAN office for more information.



ARE YOU LOOKING FOR A FAST AND AUTOMATED SOLUTION FOR SAMPLE PREPARATION? WE CAN HELP!

Automation saves time and resources. With our reagents and devices this can be done easily and reliably! By adapting the workflow, an implementation on devices such as Resolvex® A200 can be easily accomplished.



The combined use of the Steroid Panel LC-MS, process script and instrument has to be validated individually on site by each laboratory.

All pictures of the products are for illustration purpose only. Actual products may vary due to product enhancement.

Manufactured by

IBL International GmbH
Flughafenstrasse 52a
22335 Hamburg
Germany

Phone: +49 (0)40-53 28 91-0
Fax: +49 (0)40-53 28 91-11
Email: IBL@Tecan.com
www.tecan.com/immunoassays

Australia +61 3 9647 4100 **Austria** +43 62 46 89 330 **Belgium** +32 15 42 13 19 **China** +86 21 220 63 206 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170
Italy +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 17 4 **Nordic** +46 8 750 39 40 **Singapore** +65 644 41 886 **Spain** +34 93 93 595 25 31
Switzerland +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +41 44 922 81 11

Tecan Group Ltd. makes every effort to include accurate and up-to-date information; however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided. Changes can be made at any time without notice. All mentioned trademarks are protected by law. In general, the trademarks and designs referenced herein are trademarks, or registered trademarks, of Tecan Group Ltd., Maennedorf, Switzerland. A complete list may be found at www.tecan.com/trademarks. Product names and company names that are not contained in the list but are noted herein may be the trademarks of their respective owners. For technical details and detailed procedures of the specifications provided please contact your Tecan representative. This may contain reference to applications and products which are not available in all markets. Please check with your local sales representative.

© 2022 Tecan Trading AG, Switzerland, all rights reserved. For disclaimer and trademarks please visit www.tecan.com

www.tecan.com

 **TECAN.**

